




SLIT3 Antibody

Product Code	CSB-PA021769LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75094
Immunogen	Recombinant Human Slit homolog 3 protein (1199-1494AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Slit homolog 3 protein (Slit-3) (Multiple epidermal growth factor-like domains protein 5) (Multiple EGF-like domains protein 5), SLIT3, KIAA0814 MEGF5 SLIL2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLIT3

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA021769LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.