



SLIT3 Antibody, HRP conjugated

Product Code	CSB-PA021769LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75094
Immunogen	Recombinant Human Slit homolog 3 protein (1199-1494AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Slit homolog 3 protein (Slit-3) (Multiple epidermal growth factor-like domains protein 5) (Multiple EGF-like domains protein 5), SLIT3, KIAA0814 MEGF5 SLIL2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLIT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.