



SLITRK5 Antibody, HRP conjugated

Product Code	CSB-PA021774LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94991
Immunogen	Recombinant Human SLIT and NTRK-like protein 5 protein (617-890AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SLIT and NTRK-like protein 5 (Leucine-rich repeat-containing protein 11), SLITRK5, KIAA0918 LRRC11
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLITRK5