



SOCS7 Antibody, HRP conjugated

Product Code	CSB-PA022396LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14512
Immunogen	Recombinant Human Suppressor of cytokine signaling 7 protein (6-165AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Suppressor of cytokine signaling 7 (SOCS-7) (Nck, Ash and phospholipase C gamma-binding protein) (Nck-associated protein 4) (NAP-4), SOCS7, NAP4 SOCS6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SOCS7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.