



STXBP3 Antibody, HRP conjugated

Product Code	CSB-PA022907LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00186
Immunogen	Recombinant Human Syntaxin-binding protein 3 protein (266-528AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syntaxin-binding protein 3 (Platelet Sec1 protein) (PSP) (Protein unc-18 homolog 3) (Unc18-3) (Protein unc-18 homolog C) (Unc-18C), STXBP3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	STXBP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.