



TH Antibody, HRP conjugated

Product Code	CSB-PA023470LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07101
Immunogen	Recombinant Human Tyrosine 3-monooxygenase protein (2-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine 3-monooxygenase (EC 1.14.16.2) (Tyrosine 3-hydroxylase) (TH), TH, TYH
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.