





TMEM11 Antibody

Product Code	CSB-PA023679LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17152
Immunogen	Recombinant Human Transmembrane protein 11, mitochondrial protein (1-83AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 11, mitochondrial (Protein PM1) (Protein PMI), TMEM11, C17orf35 PM1
Immunogen Species	Homo sapiens (Human)
Research Area	Metabolism
Target Names	TMEM11

Image



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA023679LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lymph node tissue using CSB-PA023679LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.