



TMEM11 Antibody, Biotin conjugated

Product Code	CSB-PA023679LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17152
Immunogen	Recombinant Human Transmembrane protein 11, mitochondrial protein (1-83AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 11, mitochondrial (Protein PM1) (Protein PMI), TMEM11, C17orf35 PM1
Immunogen Species	Homo sapiens (Human)
Research Area	Metabolism
Target Names	TMEM11