



TMEM173 Antibody, Biotin conjugated

Product Code	CSB-PA023754LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WV6
Immunogen	Recombinant Human Stimulator of interferon genes protein (207-379AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Stimulator of interferon genes protein (hSTING) (Endoplasmic reticulum interferon stimulator) (ERIS) (Mediator of IRF3 activation) (hMITA) (Transmembrane protein 173), TMEM173, ERIS MITA STING
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TMEM173
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.