




TMEM43 Antibody

Product Code	CSB-PA023840LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTV4
Immunogen	Recombinant Human Transmembrane protein 43 protein (178-312AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 43 (Protein LUMA), TMEM43
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TMEM43

Image



Western Blot

Positive WB detected in: A549 whole cell lysate, Mouse brain tissue, Mouse skeletal muscle tissue


All lanes: TMEM43 antibody at 4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 45 kDa

Observed band size: 45 kDa



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA023840LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.