



TMEM43 Antibody, HRP conjugated

Product Code	CSB-PA023840LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTV4
Immunogen	Recombinant Human Transmembrane protein 43 protein (178-312AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 43 (Protein LUMA), TMEM43
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TMEM43