



TNIK Antibody, Biotin conjugated

Product Code	CSB-PA023998LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKE5
Immunogen	Recombinant Human TRAF2 and NCK-interacting protein kinase protein (523-714AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TRAF2 and NCK-interacting protein kinase (EC 2.7.11.1), TNIK, KIAA0551
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TNIK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.