



TPD52L2 Antibody, FITC conjugated

Product Code	CSB-PA024098LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43399
Immunogen	Recombinant Human Tumor protein D54 protein (1-206AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor protein D54 (hD54) (Tumor protein D52-like 2), TPD52L2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TPD52L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.