



TTK Antibody, FITC conjugated

Product Code	CSB-PA025251LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33981
Immunogen	Recombinant Human Dual specificity protein kinase TTK protein (370-492AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein kinase TTK (EC 2.7.12.1) (Phosphotyrosine picked threonine-protein kinase) (PYT), TTK, MPS1 MPS1L1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TTK