



TTK Antibody, Biotin conjugated

Product Code	CSB-PA025251LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33981
Immunogen	Recombinant Human Dual specificity protein kinase TTK protein (370-492AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein kinase TTK (EC 2.7.12.1) (Phosphotyrosine picked threonine-protein kinase) (PYT), TTK, MPS1 MPS1L1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TTK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.