



TXK Antibody, Biotin conjugated

Product Code	CSB-PA025362LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42681
Immunogen	Recombinant Human Tyrosine-protein kinase TXK protein (32-83AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase TXK (EC 2.7.10.2) (Protein-tyrosine kinase 4) (Resting lymphocyte kinase), TXK, PTK4 RLK
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TXK