



TXK Antibody, Biotin conjugated

Product Code	CSB-PA025362LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42681
Immunogen	Recombinant Human Tyrosine-protein kinase TXK protein (32-83AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase TXK (EC 2.7.10.2) (Protein-tyrosine kinase 4) (Resting lymphocyte kinase), TXK, PTK4 RLK
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TXK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.