



TYK2 Antibody, FITC conjugated

Product Code	CSB-PA025391LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29597
Immunogen	Recombinant Human Non-receptor tyrosine-protein kinase TYK2 protein (286-372AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-receptor tyrosine-protein kinase TYK2 (EC 2.7.10.2), TYK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TYK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.