



VEGFB Antibody, FITC conjugated

Product Code	CSB-PA025834LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49765
Immunogen	Recombinant Human Vascular endothelial growth factor B protein (22-207AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vascular endothelial growth factor B (VEGF-B) (VEGF-related factor) (VRF), VEGFB, VRF
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	VEGFB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.