





VPS4B Antibody

Product Code	CSB-PA025915LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75351
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 4B protein (1-117AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:200-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 4B (EC 3.6.4.6) (Cell migration-inducing gene 1 protein) (Suppressor of K(+) transport growth defect 1) (Protein SKD1), VPS4B, SKD1 VPS42
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS4B


Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA025915LA01HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA025915LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Western Blot

Positive WB detected in: Mouse heart tissue, Mouse brain tissue


All lanes: VPS4B antibody at 2.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 50 kDa

Observed band size: 50 kDa



Western Blot

Positive WB detected in: MCF-7 whole cell lysate

All lanes: VPS4B antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 50 kDa

Observed band size: 50 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.