



VPS4B Antibody, Biotin conjugated

Product Code	CSB-PA025915LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75351
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 4B protein (1-117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 4B (EC 3.6.4.6) (Cell migration-inducing gene 1 protein) (Suppressor of K(+) transport growth defect 1) (Protein SKD1), VPS4B, SKD1 VPS42
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS4B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.