




wnt8b Antibody

Product Code	CSB-PA026144LA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51029
Immunogen	Recombinant Zebrafish Protein Wnt-8b protein (24-358AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Wnt-8b, wnt8b, wnt-8b
Immunogen Species	Zebrafish
Research Area	Others
Target Names	wnt8b

Image



Western Blot

Positive WB detected in: Zebrafish tissue 40µg, 20µg, 10µg

All lanes: wnt8b antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 41 kDa

Observed band size: 41 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.