



ZNF207 Antibody, HRP conjugated

Product Code	CSB-PA026578LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43670
Immunogen	Recombinant Human BUB3-interacting and GLEBS motif-containing protein ZNF207 protein (206-463AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BUB3-interacting and GLEBS motif-containing protein ZNF207 (BuGZ) (hBuGZ) (Zinc finger protein 207), ZNF207, BUGZ
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF207
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.