



ZNF207 Antibody, Biotin conjugated

Product Code	CSB-PA026578LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43670
Immunogen	Recombinant Human BUB3-interacting and GLEBS motif-containing protein ZNF207 protein (206-463AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BUB3-interacting and GLEBS motif-containing protein ZNF207 (BuGZ) (hBuGZ) (Zinc finger protein 207), ZNF207, BUGZ
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF207