





SUMO3 Antibody

Product Code	CSB-PA05509A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55854
Immunogen	Recombinant Human Small ubiquitin-related modifier 3 protein (1-92AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Small ubiquitin-related modifier 3 (SUMO-3) (SMT3 homolog 1) (SUMO-2) (Ubiquitin-like protein SMT3A) (Smt3A), SUMO3, SMT3A SMT3H1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SUMO3

Image



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA05509A0Rb at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA05509A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.