



IL23A Antibody, Biotin conjugated

Product Code	CSB-PA06719D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPF7
Immunogen	Recombinant Human Interleukin-23 subunit alpha protein (20-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-23 subunit alpha (IL-23 subunit alpha) (IL-23-A) (Interleukin-23 subunit p19) (IL-23p19), IL23A, SGRF
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL23A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.