



MLN Antibody, FITC conjugated

Product Code	CSB-PA11109C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12872
Immunogen	Recombinant Human Promotilin protein (26-115AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Promotilin [Cleaved into: Motilin; Motilin-associated peptide (MAP)], MLN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MLN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.