



GAD2 Antibody, HRP conjugated

Product Code	CSB-PA11159B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05329
Immunogen	Recombinant Human Glutamate decarboxylase 2 protein (1-82AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate decarboxylase 2 (EC 4.1.1.15) (65 kDa glutamic acid decarboxylase) (GAD-65) (Glutamate decarboxylase 65 kDa isoform), GAD2, GAD65
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GAD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.