



GAD2 Antibody, Biotin conjugated

Product Code	CSB-PA11159D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05329
Immunogen	Recombinant Human Glutamate decarboxylase 2 protein (1-82AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate decarboxylase 2 (EC 4.1.1.15) (65 kDa glutamic acid decarboxylase) (GAD-65) (Glutamate decarboxylase 65 kDa isoform), GAD2, GAD65
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GAD2