



KRAS Antibody, FITC conjugated

Product Code	CSB-PA13967C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01116
Immunogen	Recombinant Human GTPase KRas protein (2-168AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTPase KRas (K-Ras 2) (Ki-Ras) (c-K-ras) (c-Ki-ras) [Cleaved into: GTPase KRas, N-terminally processed], KRAS, KRAS2 RASK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRAS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.