



HLA-DQB1 Antibody, Biotin conjugated

Product Code	CSB-PA14849D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01920
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DQ beta 1 chain protein (128-175AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DQ beta 1 chain (MHC class II antigen DQB1), HLA-DQB1, HLA-DQB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HLA-DQB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.