




acta1b Antibody

Product Code	CSB-PA20477A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6XNL8
Immunogen	Recombinant Zebrafish Actin protein (1-168AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	acta1b
Immunogen Species	Zebrafish
Research Area	Others
Target Names	acta1b

Image



Western Blot

Positive WB detected in: Zebrafish tissue, Zebrafish heart tissue, Zebrafish brain tissue, Zebrafish eye tissue

All lanes: acta1b antibody at 3.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 42 kDa

Observed band size: 42 kDa