



tuba8l2 Antibody, Biotin conjugated

Product Code	CSB-PA20829D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PC95
Immunogen	Recombinant Zebrafish Tubulin alpha chain protein (1-449AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin alpha chain, tuba8l2
Immunogen Species	Zebrafish
Research Area	Others
Target Names	tuba8l2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.