




sox10 Antibody

Product Code	CSB-PA20859A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q90XD1
Immunogen	Recombinant Zebrafish HMG-box transcription factor Sox10 protein (1-485AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	sox10
Immunogen Species	Zebrafish
Research Area	Others
Target Names	sox10

Image



Western Blot

Positive WB detected in: Recombinant protein
40µg, 20µg, 10µg, 5µg

All lanes: sox10 antibody at 3.4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 73 kDa

Observed band size: 73 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.