



mitfa Antibody, Biotin conjugated

Product Code	CSB-PA20909D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9PWC2
Immunogen	Recombinant Zebrafish Microphthalmia-associated transcription factor a protein (1-412AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	mitfa
Immunogen Species	Zebrafish
Research Area	Others
Target Names	mitfa
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.