



tbx16 Antibody, Biotin conjugated

Product Code	CSB-PA20979D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q504I8
Immunogen	Recombinant Zebrafish T-box gene 16 protein (1-470AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	tbx16
Immunogen Species	Zebrafish
Research Area	Others
Target Names	tbx16