



her6 Antibody, FITC conjugated

Product Code	CSB-PA20989C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PBX3
Immunogen	Recombinant Zebrafish Hairy-related 6 protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	her6
Immunogen Species	Zebrafish
Research Area	Others
Target Names	her6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.