



Larval cuticle protein A1A Antibody, Biotin conjugated

Product Code	CSB-PA304533LD01TBF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80681
Immunogen	Recombinant Tenebrio molitor Larval cuticle protein A1A protein (1-174AA)
Raised In	Rabbit
Species Reactivity	Tenebrio molitor
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Larval cuticle protein A1A (TM-LCP A1A) (TM-A1A)
Immunogen Species	Tenebrio molitor (Yellow mealworm beetle)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.