



Perlwapin Antibody, FITC conjugated

Product Code	CSB-PA306283LC01HAZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P84811
Immunogen	Recombinant <i>Haliotis laevigata</i> Perlwapin protein (1-134AA)
Raised In	Rabbit
Species Reactivity	<i>Haliotis laevigata</i>
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Perlwapin
Immunogen Species	<i>Haliotis laevigata</i>
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.