



Perlwapin Antibody, Biotin conjugated

Product Code	CSB-PA306283LD01HAZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P84811
Immunogen	Recombinant Haliotis laevigata Perlwapin protein (1-134AA)
Raised In	Rabbit
Species Reactivity	Haliotis laevigata
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Perlwapin
Immunogen Species	Haliotis laevigata
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.