




dps Antibody

Product Code	CSB-PA314354HA01MVW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C558
Immunogen	Recombinant Mycobacterium smegmatis DNA protection during starvation protein (1-183AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium smegmatis
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA protection during starvation protein (EC 1.16.-.-), dps
Immunogen Species	Mycobacterium smegmatis
Research Area	Others
Target Names	dps

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: dps antibody at 2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 36 kDa

Observed band size: 36 kDa