



N Antibody, Biotin conjugated

Product Code	CSB-PA318401LD01LNP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13699
Immunogen	Recombinant Lassa virus Nucleoprotein protein (1-569AA)
Raised In	Rabbit
Species Reactivity	Lassa virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoprotein (Nucleocapsid protein) (Protein N), N
Immunogen Species	Lassa virus
Research Area	Others
Target Names	N
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.