



LMP2 Antibody, HRP conjugated

Product Code	CSB-PA321086HB01EFA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13285
Immunogen	Recombinant Epstein-Barr virus Latent membrane protein 2 protein (1-147AA)
Raised In	Rabbit
Species Reactivity	Epstein-Barr virus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Latent membrane protein 2 (Terminal protein), LMP2
Immunogen Species	Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)
Research Area	Others
Target Names	LMP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.