




yopM Antibody

Product Code	CSB-PA324657HA01YAS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17778
Immunogen	Recombinant Yersinia pestis Outer membrane protein YopM protein (1-409AA)
Raised In	Rabbit
Species Reactivity	Yersinia pestis
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Yersinia pestis
Research Area	Others
Target Names	yopM

Image



Western Blot

Positive WB detected in: yopM 293T Transfected lysate, 293T non-Transfected lysate

All lanes: yopM antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 74 kDa

Observed band size: 74 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.