



yopM Antibody, HRP conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA324657HB01YAS |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P17778 |
| Immunogen | Recombinant Yersinia pestis Outer membrane protein YopM protein (1-409AA) |
| Raised In | Rabbit |
| Species Reactivity | Yersinia pestis |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Outer membrane protein YopM, yopM, yop48 |
| Immunogen Species | Yersinia pestis |
| Research Area | Others |
| Target Names | yopM |