



yopM Antibody, FITC conjugated

Product Code	CSB-PA324657HC01YAS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17778
Immunogen	Recombinant Yersinia pestis Outer membrane protein YopM protein (1-409AA)
Raised In	Rabbit
Species Reactivity	Yersinia pestis
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer membrane protein YopM, yopM, yop48
Immunogen Species	Yersinia pestis
Research Area	Others
Target Names	yopM