



yopM Antibody, Biotin conjugated

Product Code	CSB-PA324657HD01YAS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17778
Immunogen	Recombinant Yersinia pestis Outer membrane protein YopM protein (1-409AA)
Raised In	Rabbit
Species Reactivity	Yersinia pestis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer membrane protein YopM, yopM, yop48
Immunogen Species	Yersinia pestis
Research Area	Others
Target Names	yopM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.