




SDC4 Antibody

Product Code	CSB-PA327046HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31431
Immunogen	Recombinant Human Syndecan-4 protein (19-145AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syndecan-4 (SYND4) (Amphiglycan) (Ryudocan core protein), SDC4
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SDC4

Image



Western Blot

Positive WB detected in: Mouse heart tissue,
Mouse kidney tissue


All lanes: SDC4 antibody at 3µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 22, 17 kDa


Observed band size: 22 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA327046HA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA327046HA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA327046HA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.