



SDC4 Antibody, Biotin conjugated

Product Code	CSB-PA327046HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31431
Immunogen	Recombinant Human Syndecan-4 protein (19-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syndecan-4 (SYND4) (Amphiglycan) (Ryudocan core protein), SDC4
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SDC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.