




Cyclomaltodextrin glucanotransferase Antibody

Product Code	CSB-PA327084LA01EQJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31835
Immunogen	Recombinant Bacillus macerans Cyclomaltodextrin glucanotransferase protein (608-713AA)
Raised In	Rabbit
Species Reactivity	Bacillus macerans
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclomaltodextrin glucanotransferase (EC 2.4.1.19) (Cyclodextrin-glycosyltransferase) (CGTase)
Immunogen Species	Bacillus macerans
Research Area	Others

Image



Western Blot

Positive WB detected in Recombinant protein
All lanes: Cyclomaltodextrin glucanotransferase antibody at 2.5µg/ml
Secondary
Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 18 kDa
Observed band size: 18 kDa